



ISOTOPE
ROSATOM

Actinium-225

Form
Nitrate in solid form

Total activity
On customer's request

Specific activity
58 000 Ci/g

Half-life
10.0 days

Purity
Ac-225 > 99.9%
Ra-225 < 0.02%
Ra-224 < 0.002%
Other isotopes: < 0.0007%
Non-active impurities: < 10 mkg/mCi

Storage and transportation conditions
Ambient temperature
Type A package

Primary package
Standard glass vial for medical product

KEY ADVANTAGES:



Experienced production facility (IPPE)



Worldwide logistics



Wide fleet of certified containers

For pricing and availability please inquire

**ISOTOPE JSC: THE CORE OF
ATOMIC TECHNOLOGY
FOR LONGER LIVES**